

Anti-Trim39 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Hcm9

Item Number: CM0027402



Introduction

TRIM39 Rabbit Polyclonal Antibody validated for WB and ELISA. Detects human TRIM39, an E3 ubiquitin ligase involved in apoptosis, cell cycle regulation, and NF-κB signaling. Ideal for cancer and cell biology research.

[Learn More](#)

Parameter	Value
Product Name	TRIM39 Rabbit pAb
Remarks/Alias	TFP, RNF23, TRIM39B, TRIM39
Species	Human
GeneID (Human)	56658
GeneID	56658
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 350-450 of human TRIM39 (NP_067076.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9HCM9
Protein Weight	60kDa
Shipping	Ice bag