

Anti-Trim41 Polyclonal Antibody For Wb And Elisa - Q8Wv44

Item Number: CM0025755



Introduction

Rabbit polyclonal antibody against human TRIM41 (RINCK). Validated for WB and ELISA, cross-reacts with human, mouse, and rat. Ideal for studying innate antiviral response and ubiquitin-mediated degradation.

[Learn More](#)

Parameter	Value
Product Name	TRIM41 Rabbit pAb
Remarks/Alias	RINCK; TRIM41
Species	Human
Gene ID (Human)	90933
Gene ID	90933
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 421-630 of human TRIM41 (NP_291027.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8WV44
Protein Weight	72kDa
Shipping	Ice bag