

Anti-Trim44 Rabbit Polyclonal Antibody For Wb, Elisa - Q96Dx7

Item Number: CM0025646

Introduction

Rabbit anti-TRIM44 polyclonal antibody, validated for WB and ELISA applications, with cross-reactivity to human, mouse, and rat. Suitable for studying neuronal differentiation and PAX6 signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	TRIM44 Rabbit pAb
Remarks/Alias	AN3; MC7; DIPB; HSA249128; TRIM44
Species	Human
Gene ID (Human)	54765
Gene ID	54765
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 205-344 of human TRIM44 (NP_060053.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96DX7
Protein Weight	38kDa
Shipping	Ice bag