

Anti-Trim45 Rabbit Polyclonal Antibody For If/Icc And Elisa - Q9H8W5

Item Number: CM0018216



Introduction

Rabbit polyclonal antibody against human TRIM45 (RNF99). Validated for IF/ICC and ELISA, cross-reacts with mouse. Target plays key role in inflammatory response and ubiquitination pathways.

[Learn More](#)

Parameter	Value
Product Name	TRIM45 Rabbit pAb
Remarks/Alias	RNF99; TRIM45
Species	Human
GeneID (Human)	80263
GeneID	80263
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 100-400 of human TRIM45 (NP_001139107.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Mouse
SWISS	Q9H8W5
Protein Weight	64kDa
Shipping	Ice bag