

Anti-Trim59 Polyclonal Antibody For Wb, Elisa - Q8lwr1

Item Number: CM0013000

Introduction

Rabbit polyclonal antibody targeting human TRIM59 (TSBF1/RNF104), validated for WB and ELISA. Cross-reacts with mouse and rat. Ideal for studying ER-associated ubiquitination, innate immunity, and autophagy pathways.

[Learn More](#)

Parameter	Value
Product Name	TRIM59 Rabbit pAb
Remarks/Alias	MRF1; TSBF1; IFT80L; RNF104; TRIM57; TRIM59
Species	Human
GeneID (Human)	286827
GeneID	286827
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-330 of human TRIM59 (NP_775107.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8IWR1
Protein Weight	47kDa
Shipping	Ice bag