

Anti-Trim60 Rabbit Polyclonal Antibody For Wb And Elisa - Q495X7

Item Number: CM0021671

Introduction

Anti-TRIM60 rabbit polyclonal antibody validated for WB and ELISA. Detects TRIM60, an E3 SUMO-protein ligase regulating NF-κB signaling via TAB2 SUMOylation. Cross-reacts with mouse. 55 kDa.

[Learn More](#)

Parameter	Value
Product Name	TRIM60 Rabbit pAb
Remarks/Alias	RNF33; RNF129; TRIM60
Species	Human
GeneID (Human)	166655
GeneID	166655
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human TRIM60 (NP_689833.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q495X7
Protein Weight	55kDa
Shipping	Ice bag