

Anti-Trim62 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Bvg3

Item Number: CM0023961

Introduction

High-quality rabbit polyclonal antibody against human TRIM62 (DEAR1), an E3 ubiquitin ligase in antifungal immunity. Validated for WB, IHC-P, ELISA; cross-reacts with mouse and rat. Ideal for ubiquitin research.

[Learn More](#)

Parameter	Value
Product Name	TRIM62 Rabbit pAb
Remarks/Alias	DEAR1; TRIM62
Species	Human
Gene ID	55223
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 120-280 of human TRIM62 (NP_060677.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BVG3
Protein Weight	54kDa
Shipping	Ice bag