

Anti-Trim63 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q969Q1

Item Number: CM0026969



Introduction

TRIM63 Rabbit Polyclonal Antibody validated for WB, IHC-P, IF/ICC, and ELISA. Detects human, mouse, and rat TRIM63, a muscle-specific E3 ubiquitin ligase involved in protein degradation and atrophy.

[Learn More](#)

Parameter	Value
Product Name	TRIM63 Rabbit pAb
Remarks/Alias	IRF; SMRZ; MURF1; MURF2; RNF28; TRIM63
Species	Human
GeneID (Human)	84676
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 253-353 of human TRIM63 (NP_115977.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q969Q1
Protein Weight	40 kDa
Shipping	Ice bag