


Anti-Trim66 Rabbit Polyclonal Antibody For Wb - O15016

Item Number: CM0024176



Introduction

High-quality rabbit polyclonal antibody against human TRIM66 (TIF1D, TIF1DELTA). Suitable for WB and ELISA. Targets a nuclear transcription repressor with bromodomain and zinc-finger motifs. O15016.

[Learn More](#)

Parameter	Value
Product Name	TRIM66 Rabbit pAb
Remarks/Alias	TIF1D; C11orf29; TIF1DELTA; TRIM66
Species	Human
Gene ID (Human)	9866
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 750-970 of human TRIM66 (NP_055633.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O15016
Protein Weight	149kDa
Shipping	Ice bag