


Anti-Trim69 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q86Wt6

Item Number: CM0030420



Introduction

Rabbit polyclonal antibody against human TRIM69, validated for WB, IHC-P, IF/ICC, and ELISA. Detects endogenous TRIM69 in human and mouse samples. Suitable for ubiquitination, antiviral immunity, and centrosome dynamics research.

[Learn More](#)

Parameter	Value
Product Name	TRIM69 Rabbit pAb
Remarks/Alias	Trif; HSD34; RNF36; HSD-34; TRIM69
Species	Human
Gene ID (Human)	140691
Gene ID	140691
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-341 of human TRIM69 (NP_542783.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q86WT6
Protein Weight	57kDa
Shipping	Ice bag