

Anti-Trim72 Rabbit Polyclonal Antibody For Wb, Elisa - Q1Xh17

Item Number: CM0022674



Introduction

Rabbit polyclonal antibody targeting TRIM72 (MG53), a muscle-specific E3 ubiquitin ligase involved in membrane repair. Suitable for WB and ELISA, with cross-reactivity to human and mouse. Ideal for skeletal and cardiac muscle research.

[Learn More](#)

Parameter	Value
Product Name	TRIM72 Rabbit pAb
Remarks/Alias	MG53; TRIM72
Species	Mouse
GeneID (Human)	493829
GeneID	434246
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 200-300 of mouse TRIM72 (NP_001073401.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q1XH17
Protein Weight	53kDa
Shipping	Ice bag