

Anti-Trimethyl-Histone H3-K9 (Q16695/P68431) Rabbit Monoclonal Antibody For Epigenetics

Item Number: CM0000400



Introduction

TriMethyl-Histone H3-K9 rabbit monoclonal antibody validated for WB, ChIP-seq, CUT&Tag, DB, IF/ICC, and ELISA. Detects histone H3 trimethylated at lysine 9 across human, mouse, rat, and other species. Suited for chromatin and epigenetic studies.

[Learn More](#)

Parameter	Value
Product Name	TriMethyl-Histone H3-K9 Rabbit mAb
Remarks/Alias	H3.4; H3/g; H3FT; H3t; HIST3H3; Histone H3; HIST1H3A; TriMethyl-Histone H3-K9
Species	Human
GeneID (human)	8290/8350
GeneID	8290/8350
Immunogen	Synthetic peptide: A synthetic trimethylated peptide around K9 of human Histone H3 (NP_003520.1)
Source	Rabbit
Category	Methylated Antibodies
Application	WB, DB, IF/ICC, ELISA, ChIP, ChIP-seq, CUT&Tag
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15 kDa
Shipping	Ice bag