

Anti-Trip10 Polyclonal Antibody For Wb, If/Icc, Elisa - Q15642

Item Number: CM0025578



Introduction

High-quality rabbit polyclonal anti-TRIP10 antibody (Cdc42-interacting protein 4) validated for WB, IF/ICC, ELISA. Cross-reacts with human and mouse. Ideal for studying endocytosis, actin cytoskeleton, and GLUT4 translocation research.

[Learn More](#)

Parameter	Value
Product Name	TRIP10 Rabbit pAb
Remarks/Alias	STP; CIP4; HSTP; STOT; TRIP-10; TRIP10
Species	Human
GeneID (Human)	9322
GeneID	9322
Immunogen	Recombinant protein (Recombinant fusion protein containing a.a. 246-545 of human TRIP10)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q15642
Protein Weight	68kDa
Shipping	Ice bag