


Anti-Trip12 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q14669

Item Number: CM0029862



Introduction

Rabbit polyclonal antibody to human TRIP12 (Q14669), validated for WB, IHC-P, ELISA. Detects human, mouse, rat TRIP12. Key reagent for research on DNA repair and the ubiquitin-proteasome pathway.

[Learn More](#)

Parameter	Value
Product Name	TRIP12 Rabbit pAb
Remarks/Alias	ULF; MRD49; TRIPC; TRIP-12; TRIP12
Species	Human
GeneID (Human)	9320
GeneID	9320
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human TRIP12 (NP_001271143.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14669
Protein Weight	220kDa
Shipping	Ice bag