


Anti-Trip4 Rabbit Polyclonal Antibody For Wb, Ihc-P, Ip, Elisa - Q15650

Item Number: CM0016144



Introduction

Rabbit polyclonal antibody validated for WB, IHC-P, IP, ELISA; targets human TRIP4 (ASC-1), a transcription coactivator and RQC complex component. Cross-reacts with human, mouse, rat. Recombinant immunogen (aa 282-581).

[Learn More](#)

Parameter	Value
Product Name	TRIP4 Rabbit pAb
Remarks/Alias	ASC1; ASC-1; MDCDC; SMABF1; ZC2HC5; HsT17391; TRIP4
Species	Human
Gene ID (Human)	9325
Gene ID	9325
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 282-581 of human TRIP4 (NP_057297.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15650
Protein Weight	66kDa
Shipping	Ice bag