

Anti-Trk Fused Gene (Tfg) Rabbit Monoclonal Antibody For Wb, Elisa - Q92734

Item Number: CM0003090



Introduction

Rabbit monoclonal antibody targeting human TFG (TRK-fused gene). Validated for WB and ELISA with cross-reactivity to human and mouse. Suitable for studying ER dynamics and secretory cargo trafficking.

[Learn More](#)

Parameter	Value
Product Name	TRK fused gene (TFG) Rabbit mAb
Remarks/Alias	TF6; HMSNP; SPG57; TRKT3; TRK fused gene (TFG)
Species	Human
GeneID (Human)	10342
GeneID	10342
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q92734
Protein Weight	43kDa
Shipping	Ice bag