

Anti-Trk Fused Gene (Tfg) Rabbit Polyclonal Antibody For Wb, Elisa - Q92734

Item Number: CM0031070



Introduction

Rabbit polyclonal antibody against human TRK fused gene (TFG) protein, validated for WB and ELISA with cross-reactivity to human and mouse. Ideal for studying ER dynamics and secretory cargo trafficking.

[Learn More](#)

Parameter	Value
Product Name	TRK fused gene (TFG) Rabbit pAb
Remarks/Alias	TF6; HMSNP; SPG57; TRKT3; TRK fused gene (TFG)
Species	Human
Gene ID	10342
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human TRK fused gene (TFG) (NP_006061.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q92734
Protein Weight	43kDa
Shipping	Ice bag