


Anti-Trka+B+C Monoclonal Antibody For Wb, Ihc-P, Elisa - P04629 / Q16288 / Q16620

Item Number: CM0004031



Introduction

This recombinant rabbit monoclonal antibody recognizes TrkA, TrkB, and TrkC. Validated for Western blot, IHC-P, and ELISA applications in human, mouse, and rat tissues. Suitable for studying neurotrophin receptor signaling and neural development.

[Learn More](#)

Parameter	Value
Product Name	TrkA+B+C Rabbit mAb
Remarks/Alias	MTC; TRK; TRK1; TRKA; Trk-A; p140-TrkA; TRKC; GP145-TrkC; gp145(trkC); OBHD; TRKB; DEE58; trk-B; EIEE58; GP145-TrkB; TrkA+B+C
Species	Human
GeneID (Human)	4914/4916/4915
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 722-822 of human NTRK2 (Q16620).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P04629/Q16288/Q16620
Protein Weight	140kDa
Shipping	Ice bag