


Anti-Trkb (Gp145-Trkb) Rabbit Monoclonal Antibody For Wb And Elisa - Q16620

Item Number: CM0003678



Introduction

Rabbit monoclonal antibody targeting human TrkB (NTRK2, GP145-TrkB), validated for Western blot and ELISA. Detects mouse and rat TrkB. Ideal for neurotrophin signaling, neuronal development, and plasticity research.

[Learn More](#)

Parameter	Value
Product Name	TrkB Rabbit mAb
Remarks/Alias	OBHD; TRKB; DEE58; trk-B; EIEE58; GP145-TrkB; TrkB
Species	Human
GeneID	4915
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 32-430 of human TrkB (NP_001018074.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q16620
Protein Weight	92kDa
Shipping	Ice bag