


Anti-Trkb Rabbit Polyclonal Antibody For Wb, Elisa - Q16620

Item Number: CM0016874



Introduction

High-quality TrkB Rabbit Polyclonal Antibody for Western Blot (WB) and ELISA applications. Reacts with human and mouse TrkB (Q16620). Targets BDNF/NT-3 growth factor receptor. Reliable tool for neuroscience and neurotrophic signaling research.

[Learn More](#)

Parameter	Value
Product Name	TrkB Rabbit pAb
Remarks/Alias	OBHD; TRKB; DEE58; trk-B; EIEE58; GP145-TrkB; TrkB
Species	Human
GeneID (human)	4915
GeneID	4915
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q16620
Protein Weight	92kDa
Shipping	Ice bag