

Anti-Trkc Polyclonal Antibody For Wb, Elisa - Q16288

Item Number: CM0013351



Introduction

Rabbit polyclonal antibody targeting human TrkC (NT-3 growth factor receptor) for WB and ELISA; cross-reactivity with mouse and rat; suitable for neurotrophin signaling and neurodevelopment studies.

[Learn More](#)

Parameter	Value
Product Name	TrkC Rabbit pAb
Remarks/Alias	TRKC; GP145-TrkC; gp145(trkC); TrkC
Species	Human
Gene ID (Human)	4916
Gene ID	4916
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 40-300 of human TrkC (NP_001007157.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16288
Protein Weight	94kDa
Shipping	Ice bag