


Anti-Trkc Rabbit Monoclonal Antibody For Wb, Elisa - Q16288

Item Number: CM0006707



Introduction

High-quality rabbit monoclonal antibody targeting human TrkC (NTRK3). Validated for WB and ELISA, with cross-reactivity to mouse and rat. Suitable for neurotrophin signaling and cell differentiation research.

[Learn More](#)

Parameter	Value
Product Name	TrkC Rabbit mAb
Remarks/Alias	TRKC; GP145-TrkC; gp145(trkC); TrkC
Species	Human
Gene ID (Human)	4916
Gene ID	4916
Immunogen	Synthetic peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q16288
Protein Weight	94kDa
Shipping	Ice bag