


Anti-Trmt2A/Htf9C Polyclonal Antibody For Wb, Ihc-P, Elisa - Q8Iz69

Item Number: CM0030069



Introduction

High-quality rabbit polyclonal antibody targeting human TRMT2A/HTF9C. Validated for WB, IHC-P, ELISA, with cross-reactivity to mouse and rat. Ideal for tRNA methylation and RNA modification studies.

[Learn More](#)

Parameter	Value
Product Name	TRMT2A/HTF9C Rabbit pAb
Remarks/Alias	HTF9C; TRMT2A/HTF9C
Species	Human
Gene ID (Human)	27037
Gene ID	27037
Immunogen	Synthetic peptide - A synthetic peptide corresponding to a sequence within amino acids 500 to the C-terminus of human TRMT2A/HTF9C (NP_073564.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8Iz69
Protein Weight	69kDa
Shipping	Ice bag