


## Anti-Trmt2A/Htf9C Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q8Iz69

Item Number: CM0023012



### Introduction

High-quality rabbit polyclonal antibody against human TRMT2A/HTF9C, validated for Western blot, IHC-P, and ELISA. Targets a tRNA (uracil-5-)-methyltransferase catalyzing m5U54 formation in tRNAs and some mRNAs. Suitable for RNA modification and gene expression studies.

[Learn More](#)

| Parameter        | Value                                                                                                                   |
|------------------|-------------------------------------------------------------------------------------------------------------------------|
| Product Name     | TRMT2A/HTF9C Rabbit pAb                                                                                                 |
| Remarks/Alias    | HTF9C; TRMT2A/HTF9C                                                                                                     |
| Species          | Human                                                                                                                   |
| Gene ID          | 27037                                                                                                                   |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human TRMT2A/HTF9C (NP_073564.3) |
| Source           | Rabbit                                                                                                                  |
| Category         | Polyclonal Antibodies                                                                                                   |
| Application      | WB, IHC-P, ELISA                                                                                                        |
| Cross Reactivity | Human                                                                                                                   |
| SWISS            | Q8IZ69                                                                                                                  |
| Protein Weight   | 69kDa                                                                                                                   |
| Shipping         | Ice bag                                                                                                                 |