

Anti-Trmt2B Rabbit Pab For Wb, Ihc-P, Elisa - Q96Gj1

Item Number: CM0026584

Introduction

Rabbit polyclonal antibody against TRMT2B, a mitochondrial tRNA/rRNA methyltransferase, for WB, IHC-P, ELISA. Reacts with human and rat. Suitable for mitochondrial RNA modification and tRNA processing research.

[Learn More](#)

Parameter	Value
Product Name	TRMT2B Rabbit pAb
Remarks/Alias	CXorf34; dj341D10.3; TRMT2B
Species	Human
Gene ID (Human)	79979
Gene ID	79979
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human TRMT2B (NP_079193.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	Q96GJ1
Protein Weight	56kDa
Shipping	Ice bag