


Anti-Trmt6 Polyclonal Antibody For Wb, Elisa - Q9Uja5

Item Number: CM0028510



Introduction

Rabbit polyclonal antibody against human TRMT6, the non-catalytic subunit of tRNA m1A58 methyltransferase, also involved in mRNA N1-methylation. Validated for WB and ELISA applications.

[Learn More](#)

Parameter	Value
Product Name	TRMT6 Rabbit pAb
Remarks/Alias	GCD10; CGI-09; Gcd10p; TRMT6
Species	Human
GeneID (Human)	51605
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human TRMT6 (NP_057023.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9UJA5
Protein Weight	56kDa
Shipping	Ice bag