


Anti-Trmt61A Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - Q96Fx7

Item Number: CM0014838



Introduction

High-quality rabbit polyclonal antibody against human TRMT61A (tRNA methyltransferase catalytic subunit), validated for WB, IF/ICC, IHC-P, ELISA. Reacts with human, mouse, rat. Ideal for tRNA modification and mRNA methylation research.

[Learn More](#)

Parameter	Value
Product Name	TRMT61A Rabbit pAb
Remarks/Alias	GCD14; TRM61; Gcd14p; hTRM61; C14orf172; TRMT61A
Species	Human
Gene ID (Human)	115708
Gene ID	115708
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-289 of human TRMT61A (NP_689520.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96FX7
Protein Weight	31 kDa
Shipping	Ice bag