


Anti-Trmt61A Rabbit Monoclonal Antibody For Wb, If, Ihc, Elisa - Q96Fx7

Item Number: CM0004217



Introduction

High-quality rabbit monoclonal antibody targeting human TRMT61A. Validated for Western blot, immunofluorescence, immunohistochemistry, and ELISA. Cross-reacts with human, mouse, and rat. Suitable for tRNA modification and mRNA methylation research.

[Learn More](#)

Parameter	Value
Product Name	TRMT61A Rabbit mAb
Remarks/Alias	GCD14; TRM61; Gcd14p; hTRM61; C14orf172
Species	Human
Gene ID (Human)	115708
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-289 of human TRMT61A (NP_689520.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96FX7
Protein Weight	31kDa
Shipping	Ice bag