


Anti-Trnau1Ap Rabbit Polyclonal Antibody For Wb, Elisa - Q9Nx07

Item Number: CM0018580



Introduction

TRNAU1AP Rabbit pAb is a rabbit polyclonal antibody against human TRNAU1AP, validated for WB and ELISA. Cross-reacts with mouse. Target is a key factor in selenocysteine incorporation into selenoproteins, useful for selenoprotein synthesis and tRNA modification research.

[Learn More](#)

Parameter	Value
Product Name	TRNAU1AP Rabbit pAb
Remarks/Alias	SECP43; PRO1902; TRSPAP1; TRNAU1AP
Species	Human
Gene ID (Human)	54952
Gene ID	54952
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-287 of human TRNAU1AP (NP_060316.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9NX07
Protein Weight	32kDa
Shipping	Ice bag