


Anti-Trnt1 Polyclonal Antibody For Wb, Elisa - Q96Q11

Item Number: CM0021787



Introduction

Rabbit polyclonal antibody targeting TRNT1, a mitochondrial CCA tRNA nucleotidyltransferase. Suitable for WB and ELISA, with cross-reactivity to human, mouse, and rat. Ideal for studies in tRNA processing and mitochondrial diseases.

[Learn More](#)

Parameter	Details
Product Name	TRNT1 Rabbit pAb
Remarks/Alias	CCA1; RPEM; SIFD; MtCCA; CGI-47; TRNT1
Species	Human
Gene ID (Human)	51095
Gene ID	51095
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 135-434 of human TRNT1 (NP_886552.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96Q11
Protein Weight	50kDa
Shipping	Ice bag