


Anti-Trop-2 Monoclonal Antibody For Flow Cytometry - P09758

Item Number: CM0011032



Introduction

Rabbit anti-human TROP-2 monoclonal antibody for flow cytometry. Detects tumor-associated calcium signal transducer 2 (TACSTD2/GA733-1) in human samples. Conjugated to . Ideal for placental and pancreatic carcinoma research.

[Learn More](#)

Parameter	Specification
Product Name	Rabbit anti-Human TROP-2 mAb
Remarks/Alias	TACSTD2; EGP-1; EGP1; GA733-1; GA7331; GP50; M1S1; TROP2; tumor-associated calcium signal transducer 2
Species	Human
Gene ID (Human)	4070
Gene ID	4070
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-274 of human TROP-2 (NP_002344.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P09758
Protein Weight	35kDa
Shipping	Ice bag