

Anti-Trop-2 Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P09758

Item Number: CM0029748



Introduction

Rabbit polyclonal antibody targeting human TROP-2 (TACSTD2), validated for WB, IF/ICC, IHC-P, and ELISA. Cross-reacts with mouse and rat. Useful for tumor-associated calcium signal transduction and growth factor receptor research.

[Learn More](#)

Parameter	Value
Product Name	TROP-2 Rabbit pAb
Remarks/Alias	EGP1; GP50; M1S1; EGP-1; TROP2; GA7331; GA733-1; TROP-2
Species	Human
GeneID	4070
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-270 of human TROP-2 (NP_002344.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09758
Protein Weight	36 kDa
Shipping	Ice bag