

Anti-Trop-2 Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Mihc, Elisa - P09758

Item Number: CM0002050



Introduction

High-quality anti-TROP-2 rabbit monoclonal antibody validated for WB, IHC-P, IF-P, mIHC, and ELISA. Specific for human TROP-2 (P09758), a tumor-associated calcium signal transducer widely used in cancer research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	TROP-2 Rabbit mAb
Remarks/Alias	EGP1; GP50; M1S1; EGP-1; TROP2; GA7331; GA733-1; TROP-2
Species	Human
GeneID (Human)	4070
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-274 of human TROP-2 (NP_002344.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, mIHC, ELISA
Cross Reactivity	Human
SWISS	P09758
Protein Weight	36kDa
Shipping	Ice bag