

Anti-Trop-2 Rabbit Polyclonal Antibody For If-P And Elisa - Q8Bgv3

Item Number: CM0019850

Introduction

Rabbit polyclonal antibody against mouse TROP-2 (Tacstd2) for immunofluorescence and ELISA applications. Recognizes a cell surface glycoprotein involved in growth factor receptor signaling.

[Learn More](#)

Parameter	Value
Product Name	TROP-2 Rabbit pAb
Remarks/Alias	Ly97; EGP-1; TROP2; GA733-1
Species	Mouse
GeneID (Human)	4070
GeneID	56753
Immunogen	Recombinant protein corresponding to a sequence within amino acids 25-270 of mouse TROP-2 (NP_064431.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF-P, ELISA
Cross Reactivity	Mouse
SWISS	Q8BGV3
Protein Weight	35kDa
Shipping	Ice bag