

Anti-Tropomyosin 1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P09493

Item Number: CM0014481



Introduction

CamelBio offers rabbit polyclonal antibody against human Tropomyosin 1 (TPM1), suitable for WB, IHC-P, IF/ICC, and ELISA. Detects a 33 kDa protein; cross-reacts with mouse. Ideal for cytoskeleton and muscle contraction studies.

[Learn More](#)

Parameter	Specification
Product Name	Tropomyosin 1 Rabbit pAb
Remarks/Alias	CMH3; TMSA; CMD1Y; LVNC9; C15orf13; HEL-S-265; HTM-alpha; Tropomyosin 1
Species	Human
GeneID (Human)	7168
GeneID	7168
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-245 of human Tropomyosin 1 (NP_001018008.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P09493
Protein Weight	33kDa
Shipping	Ice bag