

Anti-Tropomyosin 1 Rabbit Monoclonal Antibody For Wb, Ip, If-P, Ihc-P, Elisa - P09493

Item Number: CM0005540



Introduction

Tropomyosin 1 Rabbit mAb for WB, IP, IF-P, IHC-P, ELISA; cross-reacts with mouse, rat. Ideal for actin-binding, muscle contraction regulation, and cytoskeleton stabilization research.

[Learn More](#)

Parameter	Value
Product Name	Tropomyosin 1 Rabbit mAb
Remarks/Alias	CMH3; TMSA; CMD1Y; LVNC9; C15orf13; HEL-S-265; HTM-alpha; Tropomyosin 1
Species	Human
Gene ID (Human)	7168
Gene ID	7168
Immunogen	Synthetic Peptide; A synthetic peptide corresponding to a sequence within amino acids 185-284 of human Tropomyosin 1 (P09493).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF-P, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P09493
Protein Weight	33 kDa
Shipping	Ice bag