

Anti-Troponin I1/Tni1 Monoclonal Antibody For Wb, Elisa - P19237

Item Number: CM0007880

Introduction

Recombinant rabbit monoclonal antibody targeting human Troponin I1 (slow skeletal muscle) for WB and ELISA. Cross-reacts with rat. Suitable for skeletal muscle research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	Troponin I1/Tni1 Rabbit mAb
Remarks/Alias	TNN1; SSTNI; Troponin I1/Tni1
Species	Human
GeneID	7135
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-187 of human Troponin I1/Tni1 (P19237)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P19237
Protein Weight	22kDa
Shipping	Ice bag