

Anti-Trove2/Ssa2 Polyclonal Antibody For Wb And Elisa - P10155

Item Number: CM0031034

Introduction

Rabbit polyclonal antibody targeting TROVE2/SSA2 (Ro60 autoantigen), validated for WB and ELISA in human, mouse, and rat. Suitable for autoimmune disease and RNA biology research.

[Learn More](#)

Parameter	Specification
Product Name	TROVE2 Rabbit pAb
Remarks/Alias	SSA2; RORNP; TROVE2
Species	Human
GeneID (Human)	6738
GeneID	6738
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human TROVE2 (NP_004591.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P10155
Protein Weight	61kDa
Shipping	Ice bag