

Anti-Troy/Tnfrsf19 Rabbit Monoclonal Antibody For Wb And Elisa - Q9Ns68

Item Number: CM0002139

Introduction

Anti-TROY/TNFRSF19 rabbit monoclonal antibody validated for Western blot and ELISA. Reacts with human, mouse, and rat TROY protein (~46 kDa). Ideal for apoptosis and NF-kappa-B signaling research.

[Learn More](#)

Parameter	Value
Product Name	TROY/TNFRSF19 Rabbit mAb
Remarks/Alias	TAJ; TROY; TRADE; TAJ-alpha; TROY/TNFRSF19
Species	Human
Gene ID (Human)	55504
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 100-200 of human TROY/TNFRSF19 (Q9NS68)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NS68
Protein Weight	46kDa
Shipping	Ice bag