


Anti-Trpc1 Polyclonal Antibody For Wb, Ihc-P, Elisa - P48995

Item Number: CM0013687



Introduction

Rabbit anti-TRPC1 polyclonal antibody validated for WB, IHC-P, and ELISA. Reacts with human and mouse TRPC1, a receptor-activated calcium channel. Ideal for ion channel research.

[Learn More](#)

Parameter	Value
Product Name	TRPC1 Rabbit pAb
Remarks/Alias	TRP1; HTRP-1; TRPC1
Species	Human
GeneID (Human)	7220
GeneID	7220
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 400-500 of human TRPC1 (NP_001238774.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P48995
Protein Weight	91kDa
Shipping	Ice bag