


Anti-Trpc1 Rabbit Monoclonal Antibody For Wb, Elisa - P48995

Item Number: CM0008349



Introduction

Rabbit monoclonal anti-TRPC1 antibody validated for Western blot and ELISA. Detects human, mouse, and rat TRPC1 (~88/91 kDa). Ideal for ion channel and calcium signaling studies.

[Learn More](#)

Parameter	Value
Product Name	TRPC1 Rabbit mAb
Remarks/Alias	TRP1; HTRP-1; TRPC1
Species	Human
Gene ID (Human)	7220
Gene ID	7220
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 400-500 of human TRPC1 (P48995)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P48995
Protein Weight	88kDa/91kDa
Shipping	Ice bag