


Anti-Trpc3 Polyclonal Antibody For Wb, Elisa - Q13507

Item Number: CM0028605



Introduction

Rabbit polyclonal antibody against human TRPC3 (TRP3, SCA41) validated for WB and ELISA. Reacts with human, mouse, rat. Target is a calcium channel involved in neuronal signaling.

[Learn More](#)

Parameter	Value
Product Name	TRPC3 Rabbit pAb
Remarks/Alias	TRP3, SCA41, TRPC3
Species	Human
GeneID	7222
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 500-600 of human TRPC3 (NP_001124170.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13507
Protein Weight	106kDa
Shipping	Ice bag