


Anti-Trpc4 Polyclonal Antibody For Wb & Elisa - Q9Ubn4

Item Number: CM0027104



Introduction

High-quality rabbit polyclonal antibody against human TRPC4. Validated for WB and ELISA, reacts with human and mouse samples. Suitable for calcium channel and endothelial cell research.

[Learn More](#)

Parameter	Value
Product Name	TRPC4 Rabbit pAb
Remarks/Alias	TRP4; HTRP4; HTRP-4; TRPC4
Species	Human
GeneID (human)	7223
GeneID	7223
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 768-977 of human TRPC4 (NP_057263.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9UBN4
Protein Weight	112kDa
Shipping	Ice bag