


## Anti-Trpc5 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9UI62

Item Number: CM0013734



### Introduction

Affinity-purified rabbit polyclonal antibody against human TRPC5, validated for Western blot (WB), immunofluorescence (IF/ICC), and ELISA. Cross-reacts with mouse. Ideal for studying calcium channels and neuronal function.

[Learn More](#)

| Parameter        | Value                                                                                                 |
|------------------|-------------------------------------------------------------------------------------------------------|
| Product Name     | TRPC5 Rabbit pAb                                                                                      |
| Remarks/Alias    | TRP5; PPP1R159; TRPC5                                                                                 |
| Species          | Human                                                                                                 |
| GeneID (Human)   | 7224                                                                                                  |
| GeneID           | 7224                                                                                                  |
| Immunogen        | Synthetic peptide corresponding to a sequence within amino acids 800-900 of human TRPC5 (NP_036603.1) |
| Source           | Rabbit                                                                                                |
| Category         | Polyclonal Antibodies                                                                                 |
| Application      | WB, IF/ICC, ELISA                                                                                     |
| Cross Reactivity | Human, Mouse                                                                                          |
| SWISS            | Q9UL62                                                                                                |
| Protein Weight   | 111kDa                                                                                                |
| Shipping         | Ice bag                                                                                               |