


Anti-Trpc6 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Y210

Item Number: CM0030032



Introduction

Rabbit polyclonal antibody against human TRPC6, validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat. Targets a calcium-permeable cation channel crucial for Gq-coupled signaling and podocyte homeostasis.

[Learn More](#)

Parameter	Specification
Product Name	TRPC6 Rabbit pAb
Remarks/Alias	TRP6; FSGS2; TRPC6
Species	Human
GeneID (Human)	7225
GeneID	7225
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 466-565 of human TRPC6 (NP_004612.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y210
Protein Weight	106 kDa
Shipping	Ice bag