

Anti-Trpm1 Rabbit Polyclonal Antibody For Wb And Elisa - Q7Z4N2

Item Number: CM0023075

Introduction

TRPM1 rabbit polyclonal antibody validated for Western Blot and ELISA. Reacts with human TRPM1, a cation channel involved in retinal photoresponse and potential metastasis suppression. Ideal for vision and oncology research.

[Learn More](#)

Parameter	Value
Product Name	TRPM1 Rabbit pAb
Remarks/Alias	MLSN1; CSNB1C; LTRPC1; TRPM1
Species	Human
Gene ID	4308
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1460-1540 of human TRPM1 (NP_002411.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q7Z4N2
Protein Weight	182kDa
Shipping	Ice bag