


Anti-Trpm2 Rabbit Polyclonal Antibody For Wb, If-P, Elisa - 094759

Item Number: CM0025537



Introduction

High-quality rabbit anti-TRPM2 polyclonal antibody validated for western blotting (WB), immunofluorescence (IF-P), and ELISA. Cross-reacts with mouse and rat TRPM2. Suitable for studying ion channels, oxidative stress, and calcium signaling.

[Learn More](#)

Parameter	Specification
Product Name	TRPM2 Rabbit pAb
Remarks/Alias	KNP3; EREG1; TRPC7; LTRPC2; NUDT9H; LTRpC-2; NUDT9L1; TRPM2
Species	Human
GeneID (Human)	7226
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1204-1503 of human TRPM2 (NP_003298.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	094759
Protein Weight	171kDa
Shipping	Ice bag