


Anti-Trpm4 Monoclonal Antibody For Wb And Elisa - Q8Td43

Item Number: CM0008700



Introduction

Rabbit monoclonal antibody against human TRPM4, validated for Western blot and ELISA. Suitable for studying calcium-activated cation channel function in cardiovascular, neuronal, and immune research.

[Learn More](#)

Parameter	Value
Product Name	TRPM4 Rabbit mAb
Remarks/Alias	EKVP6; LTrpC4; PFHB1B; TRPM4B; hTRPM4
Species	Human
Gene ID (Human)	54795
Gene ID	54795
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1097-1214 of human TRPM4 (NP_060106.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q8TD43
Protein Weight	134kDa
Shipping	Ice bag