


Anti-Trpm4 Polyclonal Antibody For Wb, Elisa - Q8Td43

Item Number: CM0015822



Introduction

TRPM4 rabbit polyclonal antibody validated for Western blot and ELISA. Recognizes human, mouse, and rat TRPM4, a calcium-activated cation channel involved in cardiac and immune function. Suitable for ion channel research.

[Learn More](#)

Parameter	Value
Product Name	TRPM4 Rabbit pAb
Remarks/Alias	EKVP6; LTrpC4; PFHB1B; TRPM4B; hTRPM4; TRPM4
Species	Human
Gene ID	54795
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human TRPM4 (NP_060106.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8TD43
Protein Weight	134kDa
Shipping	Ice bag