


Anti-Trpm8 Rabbit Monoclonal Antibody For Wb, Elisa - Q7Z2W7

Item Number: CM0008943



Introduction

Rabbit monoclonal antibody against human TRPM8 (Q7Z2W7). Validated for WB and ELISA. Suitable for cold-sensing ion channel research in prostate cancer and other tumors.

[Learn More](#)

Parameter	Value
Product Name	TRPM8 Rabbit mAb
Remarks/Alias	TRPP8; LTRPC6; trp-p8; LTrpC-6; TRPM8
Species	Human
Gene ID (Human)	79054
Gene ID	79054
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 900-1000 of human TRPM8 (Q7Z2W7).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q7Z2W7
Protein Weight	128kDa
Shipping	Ice bag